

The Golden Rules of Thermofocus

1. In order to take a body temperature, position and hold the Thermofocus *exactly in the middle of the forehead and always use the “face” button only.*
2. The Thermofocus needs to stabilize at the same ambient room temperature as the baby (or other patient) before an accurate temperature reading is possible. *You will need to wait several minutes for the Thermofocus to recalibrate if*
 - the patient has just come from a warmer or colder environment than where the Thermofocus has been stored.
 - The Thermofocus has been stored in a different room or environment than the patient.
3. An accurate reading may not be possible on a patient with a perspiring forehead, even if the forehead is dried off first.

For details on how to measure the temperature of other objects, and complete details on usage of Thermofocus, please refer to the enclosed instruction booklet.